

EBP08001 - IMMUNIZING PEPTIDE

SKU: EBP08001



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001091894.1; NP_004078.2; NP_001191316.1; NP_001191317.1; NP_001353132.1; NP_001353133.1; NP_001353140.1; NP_001353145.1; NP_001353150.1; NP_001353151.1
Immunogen	Peptide with sequence C-NSKTRDKGEIPDD, from the internal region of the protein sequence according to NP_001091894.1; NP_004078.2; NP_001191316.1; NP_001191317.1; NP_001353132.1; NP_001353133.1; NP_001353140.1; NP_001353145.1; NP_001353150.1; NP_001353151.1.
Product Comments	The peptide was used in the production of EB08001 - Goat Anti-DLG1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-NSKTRDKGEIPDD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)