

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **EBP08001 - IMMUNIZING PEPTIDE**

**SKU:** EBP08001

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Shipped at ambient temperature, store at -20°C.

NP\_001091894.1; NP\_004078.2; NP\_001191316.1; NP\_001191317.1; NP\_001353132.1; NP\_001353133.1; Accession

NP\_001353140.1; NP\_001353145.1; NP\_001353150.1; NP\_001353151.1 ID

Peptide with sequence C-NSKTRDKGEIPDD, from the internal region of the protein sequence according to

Immunogen NP 001091894.1; NP 004078.2; NP 001191316.1; NP 001191317.1; NP 001353132.1; NP 001353133.1;

NP 001353140.1; NP 001353145.1; NP 001353150.1; NP 001353151.1.

The peptide was used in the production of EB08001 - Goat Anti-DLG1 Antibody. It is currently untested in **Product** 

Comments blocking applications.

**Peptide** 

C-NSKTRDKGEIPDD Sequence

Shipping

Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade