

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP08009 - IMMUNIZING PEPTIDE**

**SKU:** EBP08009

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_571154.1 ID

Immunogen Peptide with sequence KRKIKQTVTSLPHQ-C, from the N Terminus of the protein sequence according to NP\_571154.1.

The peptide was used in the production of EB08009 - Goat Anti-Her6 (zebrafish) Antibody. It is currently untested **Product** 

**Comments** in blocking applications.

Peptide

KRKIKQTVTSLPHQ-C Sequence

Shipping

Ambient Instructions

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade