

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08025 - IMMUNIZING PEPTIDE

SKU: EBP08025

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Shipped at ambient temperature, store at -20°C.

 $NP_955392.1; \ NP_955393.1; \ NP_955394.1; \ NP_955395.1; \ NP_003279.2; \ NP_955391.1; \ NP_001230820.1; \ NP_0012308$ Accession

NP_001230821.1; NP_001230823.1 ID

Peptide with sequence C-SGDKPLSDPAP, from the C Terminus of the protein sequence according to

Immunogen NP 955392.1; NP 955393.1; NP 955394.1; NP 955395.1; NP 003279.2; NP 955391.1; NP 001230820.1;

NP_001230821.1; NP_001230823.1.

The peptide was used in the production of EB08025 - Goat Anti-TPD52L2 / D54 Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-SGDKPLSDPAP Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade