

EBP08025 - IMMUNIZING PEPTIDE

SKU: EBP08025



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_955392.1; NP_955393.1; NP_955394.1; NP_955395.1; NP_003279.2; NP_955391.1; NP_001230820.1; NP_001230821.1; NP_001230823.1
Immunogen	Peptide with sequence C-SGDKPLSDPAP, from the C Terminus of the protein sequence according to NP_955392.1; NP_955393.1; NP_955394.1; NP_955395.1; NP_003279.2; NP_955391.1; NP_001230820.1; NP_001230821.1; NP_001230823.1.
Product Comments	The peptide was used in the production of EB08025 - Goat Anti-TPD52L2 / D54 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-SGDKPLSDPAP
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)