

EBP08026 - IMMUNIZING PEPTIDE

SKU: EBP08026



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_003278.1; NP_001003395.1 ; NP_001003396.1; NP_001003397.1
Immunogen	Peptide with sequence C-EPLQGTDEDAVASAD, from the internal region of the protein sequence according to NP_003278.1; NP_001003395.1 ; NP_001003396.1; NP_001003397.1.
Product Comments	The peptide was used in the production of EB08026 - Goat Anti-TPD52L1 / D53 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-EPLQGTDEDAVASAD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)