

EBP08027 - IMMUNIZING PEPTIDE

SKU: EBP08027



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_277051.3; NP_001001874.1; NP_001001875.1
Immunogen	Peptide with sequence PHARTETSVGTYESH-C, from the N Terminus of the protein sequence according to NP_277051.3; NP_001001874.1; NP_001001875.1.
Product Comments	The peptide was used in the production of EB08027 - Goat Anti-TPD52L3 / D55 Antibody. It is currently untested in blocking applications.
Peptide Sequence	PHARTETSVGTYESH-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)