

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP08030 - IMMUNIZING PEPTIDE**

**SKU:** EBP08030

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_037481.1 ID

Peptide with sequence C-NEDGSPNLDFKP, from the internal region (near the C Terminus) of the protein sequence

according to NP\_037481.1.

The peptide was used in the production of EB08030 - Goat Anti-Neudesin / NENF Antibody. It is currently untested Product

**Comments** in blocking applications.

Peptide

C-NEDGSPNLDFKP Sequence

Shipping

Ambient Instructions

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade