

EBP08090 - IMMUNIZING PEPTIDE

SKU: EBP08090



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001035959.1; NP_001035960.1; NP_001035961.1; NP_001035962.1; NP_005126.1; NP_598408.1
Immunogen	Peptide with sequence SRHIDVCSKTKE, from the internal region of the protein sequence according to NP_001035959.1; NP_001035960.1; NP_001035961.1; NP_001035962.1; NP_005126.1; NP_598408.1.
Product Comments	The peptide was used in the production of EB08090 - Goat Anti-KCC3 / SLC12A6 Antibody. It is currently untested in blocking applications.
Peptide Sequence	SRHIDVCSKTKE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)