

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08108 - IMMUNIZING PEPTIDE

SKU: EBP08108

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

ID

Shipped at ambient temperature, store at -20°C. Instructions

Accession

 $NP_001287914.1; \ NP_001287915.1; \ NP_001287916.1; \ NP_001287917.1; \ NP_001287918.1; \ NP_067013.1; \ NP_0$

Peptide with sequence C-SKKFKGEDKMD, from the internal region of the protein sequence according to Immunogen

NP_001287914.1; NP_001287915.1; NP_001287916.1; NP_001287917.1; NP_001287918.1; NP_067013.1;.

The peptide was used in the production of EB08108 - Goat Anti-PTBP2 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-SKKFKGEDKMD Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade