

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08116 - IMMUNIZING PEPTIDE

SKU: EBP08116

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_078928.3 ID

Immunogen Peptide with sequence C-ELPTEPTDDNP, from the internal region of the protein sequence according to NP_078928.3.

The peptide was used in the production of EB08116 - Goat Anti-Leucine-rich repeat kinase 1 Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-ELPTEPTDDNP Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade