

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08135 - IMMUNIZING PEPTIDE

SKU: EBP08135

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Accession

Instructions Shipped at ambient temperature, store at -20°C.

ID NP 001337232.1; NP 001337233.1; NP 001337236.1; NP 001337243.1; NP 001337249.1

Peptide with sequence C-RREVRSPTKAPH, from the internal region of the protein sequence according to

NP 001422.1; NP 001129026.1; NP 001186317.1; NP 001239589.1; NP 001337228.1; NP 001337230.1;

Immunogen NP_001422.1; NP_001129026.1; NP_001186317.1; NP_001239589.1; NP_001337228.1; NP_001337230.1;

 $NP_001337232.1;\ NP_001337233.1;\ NP_001337236.1;\ NP_001337243.1;\ NP_001337249.1.$

Product The peptide was used in the production of EB08135 - Goat Anti-EPB41L2 / 4.1G (aa 593 to 604) Antibody. It is

Comments currently untested in blocking applications.

Peptide

Sequence

C-RREVRSPTKAPH

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade