

EBP08135 - IMMUNIZING PEPTIDE

SKU: EBP08135



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001422.1; NP_001129026.1; NP_001186317.1; NP_001239589.1; NP_001337228.1; NP_001337230.1; NP_001337232.1; NP_001337233.1; NP_001337236.1; NP_001337243.1; NP_001337249.1
Immunogen	Peptide with sequence C-RREVRSPKAPH, from the internal region of the protein sequence according to NP_001422.1; NP_001129026.1; NP_001186317.1; NP_001239589.1; NP_001337228.1; NP_001337230.1; NP_001337232.1; NP_001337233.1; NP_001337236.1; NP_001337243.1; NP_001337249.1.
Product Comments	The peptide was used in the production of EB08135 - Goat Anti-EPB41L2 / 4.1G (aa 593 to 604) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-RREVRSPKAPH
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)