

EBP08138 - IMMUNIZING PEPTIDE

SKU: EBP08138



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_001077395.1; NP_079356.3
Immunogen	Peptide with sequence EKHEKGAQKTDCQ, from the internal region of the protein sequence according to NP_001077395.1; NP_079356.3.
Product	The peptide was used in the production of EB08138 - Goat Anti-ASRGL1 / ALP Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	EKHEKGAQKTDCQ
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)