

EBP08151 - IMMUNIZING PEPTIDE

SKU: EBP08151



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_060697.4
Immunogen	Peptide with sequence C-ETTEVRKTKIRN, from the internal region of the protein sequence according to NP_060697.4.
Product Comments	The peptide was used in the production of EB08151 - Goat Anti-UBE1L2 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-ETTEVRKTKIRN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png