

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08152 - IMMUNIZING PEPTIDE

SKU: EBP08152

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

ID

 $NP_001035982.1; NP_112194.2; NP_001245295.1; NP_001245296.1; NP_001245297.1; NP_001245298.1; NP_001245299.1; NP_001245299.1;$

Peptide with sequence C-RTPMPKDVRQS, from the internal region (near the C terminus) of the protein sequence

Immunogen according to NP_001035982.1; NP_112194.2; NP_001245295.1; NP_001245296.1; NP_001245297.1;

NP_001245298.1.

The peptide was used in the production of EB08152 - Goat Anti-DIAPH3 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

Sequence

C-RTPMPKDVRQS

Shipping

Ambient Instructions

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png