

EBP08152 - IMMUNIZING PEPTIDE

SKU: EBP08152



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001035982.1; NP_112194.2; NP_001245295.1;NP_001245296.1; NP_001245297.1; NP_001245298.1 Peptide with sequence C-RTPMPKDVRQS, from the internal region (near the C terminus) of the protein sequence according to NP_001035982.1; NP_112194.2; NP_001245295.1;NP_001245296.1; NP_001245297.1; NP_001245298.1.
Immunogen	
Product Comments	The peptide was used in the production of EBP08152 - Goat Anti-DIAPH3 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-RTPMPKDVRQS
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)