

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP08178 - IMMUNIZING PEPTIDE**

**SKU:** EBP08178

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_003873.1; NP\_543028.1 ID

Immunogen Peptide with sequence C-ESYPDFSEIAN, from the C Terminus of the protein sequence according to NP\_003873.1; NP\_543028.1.

The peptide was used in the production of EB08178 - Goat Anti-WISP1 Antibody. It is currently untested in blocking **Product** 

Comments applications.

**Peptide** 

C-ESYPDFSEIAN Sequence

**Shipping** 

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_plus\_medium.png$ Grade