

EBP08199 - IMMUNIZING PEPTIDE

SKU: EBP08199



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_596867.1; NP_002202.2; NP_389647.1; NP_391987.1; NP_391989.1; NP_391988.1
Immunogen	Peptide with sequence NKCPKKDSDSKIRP, from the Internal region of the protein sequence according to NP_596867.1; NP_002202.2; NP_389647.1; NP_391987.1; NP_391989.1; NP_391988.1.
Product Comments	The peptide was used in the production of EBP08199 - Goat Anti-ITGB1 / CD29 Antibody. It is currently untested in blocking applications.
Peptide Sequence	NKCPKKDSDSKIRP
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)