

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08202 - IMMUNIZING PEPTIDE

SKU: EBP08202

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_002554.1 ID

Immunogen Peptide with sequence C-RSEANLQSKSED, from the C terminus of the protein sequence according to NP_002554.1.

Product The peptide was used in the production of EB08202 - Goat Anti-P2RY1 (aa345-356) Antibody. It is currently

Comments untested in blocking applications.

Peptide

C-RSEANLQSKSED Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade