

EBP08202 - IMMUNIZING PEPTIDE

SKU: EBP08202

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.**Unit Size** 100 µg**Storage Instructions** Shipped at ambient temperature, store at -20°C.**Accession ID** NP_002554.1**Immunogen** Peptide with sequence C-RSEANLQSKSED, from the C terminus of the protein sequence according to NP_002554.1.**Product Comments** The peptide was used in the production of EB08202 - Goat Anti-P2RY1 (aa345-356) Antibody. It is currently untested in blocking applications.**Peptide Sequence** C-RSEANLQSKSED**Shipping Instructions** Ambient**Product Grade** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)