

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08237 - IMMUNIZING PEPTIDE

SKU: EBP08237

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

 $NP_065960.2; \ NP_001171866.1; \ NP_001171868.1; \ NP_001317239.1$ ID

Peptide with sequence C-ENLPQSPGTDQHD, from the internal region of the protein sequence according to

Immunogen NP_065960.2; NP_001171866.1; NP_001171868.1; NP_001317239.1.

The peptide was used in the production of EB08237 - Goat Anti-EBPL41L5 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-ENLPQSPGTDQHD Sequence

Shipping

Instructions Ambient

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade