

EBP08237 - IMMUNIZING PEPTIDE

SKU: EBP08237



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_065960.2; NP_001171866.1; NP_001171868.1; NP_001317239.1
Immunogen	Peptide with sequence C-ENLPQSPGTDQHD, from the internal region of the protein sequence according to NP_065960.2; NP_001171866.1; NP_001171868.1; NP_001317239.1.
Product Comments	The peptide was used in the production of EB08237 - Goat Anti-EBPL41L5 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-ENLPQSPGTDQHD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)