

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP08263 - IMMUNIZING PEPTIDE**

**SKU:** EBP08263

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_004050.3; NP\_001116144.1; NP\_001274319.1; NP\_001339926.1; ID

Peptide with sequence C-DISDFKRKMPD, from the internal region of the protein sequence according to Immunogen

NP\_004050.3; NP\_001116144.1; NP\_001274319.1; NP\_001339926.1;.

The peptide was used in the production of EB08263 - Goat Anti-CCBL1 Antibody. It is currently untested in Product

Comments blocking applications.

Peptide

C-DISDFKRKMPD Sequence

Shipping

Ambient Instructions

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade