

EBP08275 - IMMUNIZING PEPTIDE

SKU: EBP08275



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_006536.3
Immunogen	Peptide with sequence QYSNVYCPTPKV, from the internal region of the protein sequence according to NP_006536.3.
Product Comments	The peptide was used in the production of EB08275 - Goat Anti-TUSC4 / NPRL2 (aa 140 to 151) Antibody. It is currently untested in blocking applications.
Peptide Sequence	QYSNVYCPTPKV
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)