

EBP08318 - IMMUNIZING PEPTIDE

SKU: EBP08318



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_937756.1; NP_937757.1; NP_443132.2
Immunogen	Peptide with sequence KQRTAPKTQYTGC, from the internal region of the protein sequence according to NP_937756.1; NP_937757.1; NP_443132.2.
Product Comments	The peptide was used in the production of EB08318 - Goat Anti-CSMD3 Antibody. It is currently untested in blocking applications.
Peptide Sequence	KQRTAPKTQYTGC
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)