

EBP08335 - IMMUNIZING PEPTIDE

SKU: EBP08335



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_079535.3 ; NP_751897.1; NP_002432.1
Immunogen	Peptide with sequence GDLHCEIRDQET, from the internal region of the protein sequence according to NP_079535.3 ; NP_751897.1; NP_002432.1.
Product Comments	The peptide was used in the production of EB08335 - Goat Anti-MSH5 Antibody. It is currently untested in blocking applications.
Peptide Sequence	GDLHCEIRDQET
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)