

EBP08340 - IMMUNIZING PEPTIDE

SKU: EBP08340



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_064524.3
Immunogen	Peptide with sequence DCEEKDWDFSMN, from the internal region of the protein sequence according to NP_064524.3.
Product Comments	The peptide was used in the production of EB08340 - Goat Anti-BDH2 / DHRS6 (aa 92 to 103) Antibody. It is currently untested in blocking applications.
Peptide Sequence	DCEEKDWDFSMN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)