

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08341 - IMMUNIZING PEPTIDE

SKU: EBP08341

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

100 μg **Unit Size**

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 065119.1

ID

Peptide with sequence NDCFSLNYPGNPCP, from the N terminus of the protein sequence according to Immunogen

NP 065119.1.

The peptide was used in the production of EB08341 - Goat Anti-HRASLS / H-REV107 (N term) Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

NDCFSLNYPGNPCP Sequence

Shipping Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png