

EBP08354 - IMMUNIZING PEPTIDE

SKU: EBP08354



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_000152.1 ; NP_001019195.1 ; NP_001019241.1 ; NP_001019242.1
Immunogen	Peptide with sequence C-GKVHIGYLPNKQ, from the internal region of the protein sequence according to NP_000152.1 ; NP_001019195.1 ; NP_001019241.1 ; NP_001019242.1 .
Product	The peptide was used in the production of EB08354 - Goat Anti-GCH1 Antibody. It is currently untested in
Comments	blocking applications.
Peptide Sequence	C-GKVHIGYLPNKQ
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)