

EBP08374 - IMMUNIZING PEPTIDE

SKU: EBP08374



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_000910.2; NP_620121.1; NP_620176.1; NP_620177.1
Immunogen	Peptide with sequence C-DQEKEDDGSESEE, from the C Terminus of the protein sequence according to NP_000910.2; NP_620121.1; NP_620176.1; NP_620177.1.
Product	The peptide was used in the production of EB08374 - Goat Anti-PAM / PAL Antibody. It is currently untested in
Comments	blocking applications.
Peptide Sequence	C-DQEKEDDGSESEE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)