

EBP08378 - IMMUNIZING PEPTIDE

SKU: EBP08378



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_878919.2; NP_001008410.1; NP_060704.2

Immunogen Peptide with sequence C-QPYVQESQNKFFK, from the internal region of the protein sequence according to NP_878919.2; NP_001008410.1; NP_060704.2.

Product Comments The peptide was used in the production of EBP08378 - Goat Anti-Dudulin 2 / STEAP3 Antibody. It is currently untested in blocking applications.

Peptide Sequence C-QPYVQESQNKFFK

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png