

EBP08380 - IMMUNIZING PEPTIDE

SKU: EBP08380



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_006756.2; NP_839952.1
Immunogen	Peptide with sequence PQRQCSVCRQAN, from the internal region of the protein sequence according to NP_006756.2; NP_839952.1.
Product	The peptide was used in the production of EB08380 - Goat Anti-TUSC3 / N33 Antibody. It is currently untested
Comments	in blocking applications.
Peptide Sequence	PQRQCSVCRQAN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)