

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP08445 - IMMUNIZING PEPTIDE**

**SKU:** EBP08445

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP\_001070671.1; NP\_065705.3; NP\_001269730.1; NP\_001269731.1 ID

The peptide was used in the production of EB08445 - Goat Anti-SENP7 (internal region) Antibody. It is currently Product

**Comments** untested in blocking applications.

Peptide

C-TRKENNLTEDNPN Sequence

Shipping

Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade