

# EBP08466 - IMMUNIZING PEPTIDE

**SKU:** EBP08466



---

## SPECIFICATIONS

**Formulation** 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 µg

**Storage Instructions** Shipped at ambient temperature, store at -20°C.

**Accession ID** NP\_006142.1

**Immunogen** Peptide with sequence C-, from the internal region of the protein sequence according to NP\_006142.1.

**Product** The peptide was used in the production of EB08466 - Goat Anti-Lactoperoxidase Antibody. It is currently

**Comments** untested in blocking applications.

**Peptide Sequence** C-

**Shipping Instructions** Ambient

**Product Grade** [https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\\_medium.png](https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png)

## DOCUMENTS

- [Data Sheet](#)