

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP08466 - IMMUNIZING PEPTIDE**

**SKU:** EBP08466

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μg

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 006142.1

Immunogen Peptide with sequence C-, from the internal region of the protein sequence according to NP\_006142.1.

**Product** The peptide was used in the production of EB08466 - Goat Anti-Lactoperoxidase Antibody. It is currently

**Comments** untested in blocking applications.

**Peptide** 

Sequence

Shipping Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade