

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08472 - IMMUNIZING PEPTIDE

SKU: EBP08472

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_009206.1 ID

Immunogen Peptide with sequence C-TKNGQKRAKLRRVH, from the interal region of the protein sequence according to NP_009206.1.

The peptide was used in the production of EB08472 - Goat Anti-TUSC2 / FUS1 Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

C-TKNGQKRAKLRRVH Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade