

EBP08493 - IMMUNIZING PEPTIDE

SKU: EBP08493



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_647536.1

Immunogen Peptide with sequence QKYCIYQHFQFQ, from the N Terminus of the protein sequence according to NP_647536.1.

Product Comments The peptide was used in the production of EB08493 - Goat Anti-CHRFAM7A / CHRNA7-FAM7A Antibody. It is currently untested in blocking applications.

Peptide Sequence QKYCIYQHFQFQ

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)