

EBP08496 - IMMUNIZING PEPTIDE

SKU: EBP08496



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_005798.3; NP_001121180.2; NP_001290161.1
Immunogen	Peptide with sequence DYESFCAEVHNP, from the internal region of the protein sequence according to NP_005798.3; NP_001121180.2; NP_001290161.1.
Product Comments	The peptide was used in the production of EBP08496 - Goat Anti-Lubricin / PRG4 Antibody. It is currently untested in blocking applications.
Peptide Sequence	DYESFCAEVHNP
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)