

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08499 - IMMUNIZING PEPTIDE

SKU: EBP08499

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_055854.1; NP_733820.1 ID

Immunogen Peptide with sequence C-RDAADTDRIMNHSS, from the internal region of the protein sequence according to NP_055854.1; NP_733820.1.

The peptide was used in the production of EB08499 - Goat Anti-NMNAT2 Antibody. It is currently untested in **Product**

Comments blocking applications. Peptide

C-RDAADTDRIMNHSS Sequence

Ambient Instructions Shipping

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade