

EBP08509 - IMMUNIZING PEPTIDE

SKU: EBP08509



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_000028.3; NP_065208.2; NP_065209.2; NP_065210.2; NP_065211.2; NP_065213.2
Immunogen	Peptide with sequence C-QIVKRASLKRKQ, from the C Terminus of the protein sequence according to NP_000028.3; NP_065208.2; NP_065209.2; NP_065210.2; NP_065211.2; NP_065213.2.
Product Comments	The peptide was used in the production of EBP08509 - Goat Anti-Ankyrin 1 / ANK1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-QIVKRASLKRKQ
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)