

# EBP08550 - IMMUNIZING PEPTIDE

**SKU:** EBP08550



---

## SPECIFICATIONS

<b>Formulation</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Unit Size</b>	100 µg
<b>Storage</b>	Shipped at ambient temperature, store at -20°C.
<b>Instructions</b>	
<b>Accession ID</b>	NP_001092692.1; NP_001092693.1; NP_079379.2; AAK77221.1
<b>Immunogen</b>	Peptide with sequence CKTIVDALHSTSGN, from the C or N Terminus of the protein sequence according to NP_001092692.1; NP_001092693.1; NP_079379.2; AAK77221.1.
<b>Product Comments</b>	The peptide was used in the production of EBP08550 - Goat Anti-CMG1 / CCDC2 / IFT74 Antibody. It is currently untested in blocking applications.
<b>Peptide Sequence</b>	CKTIVDALHSTSGN
<b>Shipping Instructions</b>	Ambient
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>

## DOCUMENTS

- [Data Sheet](#)