

EBP08552 - IMMUNIZING PEPTIDE

SKU: EBP08552



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_009213.1; NP_899237.1; NP_899239.2; NP_899240.2
Immunogen	Peptide with sequence HKFKKGDEYDVC, from the internal region of the protein sequence according to NP_009213.1; NP_899237.1; NP_899239.2; NP_899240.2.
Product	The peptide was used in the production of EB08552 - Goat Anti-RNF13 / RZF Antibody. It is currently untested
Comments	in blocking applications.
Peptide Sequence	HKFKKGDEYDVC
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)