

EBP08564 - IMMUNIZING PEPTIDE

SKU: EBP08564



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_036066.2
Immunogen	Peptide with sequence QQKPD CRPYWEKDD, from the Internal region of the protein sequence according to NP_036066.2.
Product	The peptide was used in the production of EB08564 - Goat Anti-FTO (Mouse) Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	QQKPD CRPYWEKDD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)