

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08570 - IMMUNIZING PEPTIDE

SKU: EBP08570

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession $NP_005661.1; \ NP_001018051.1; \ NP_001346986.1; \ NP_001355059.1$ ID

Peptide with sequence C-EATGHTNEMKHTTD, from the internal region of the protein sequence according to

Immunogen NP_005661.1; NP_001018051.1; NP_001346986.1; NP_001355059.1.

The peptide was used in the production of EB08570 - Goat Anti-Laforin (isoform a) Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

C-EATGHTNEMKHTTD Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade