

EBP08594 - IMMUNIZING PEPTIDE

SKU: EBP08594



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_002471.1, NP_001231919.1, NP_001231921 .1, NP_001137358.1
Immunogen	Peptide with sequence C-REDEYKQKYSRTYD, from the C Terminus of the protein sequence according to NP_002471.1, NP_001231919.1, NP_001231921 .1, NP_001137358.1.
Product Comments	The peptide was used in the production of EB08594 - Goat Anti-PPP1R12A Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-REDEYKQKYSRTYD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)