

EBP08624 - IMMUNIZING PEPTIDE

SKU: EBP08624



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001094077.1; NP_057270.1
Immunogen	Peptide with sequence QSPETPAPQNQTSR-C, from the N Terminus (near) of the protein sequence according to NP_001094077.1; NP_057270.1.
Product Comments	The peptide was used in the production of EBP08624 - Goat Anti-SERPINA10 / ZPI Antibody. It is currently untested in blocking applications.
Peptide Sequence	QSPETPAPQNQTSR-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)