

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP08658 - IMMUNIZING PEPTIDE

SKU: EBP08658

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_006668.1 ID

 $\begin{tabular}{ll} \textbf{Immunogen} & Peptide with sequence RSRNPGKDDCAKE, from the internal region of the protein sequence according to NP_006668.1. \end{tabular}$

The peptide was used in the production of EB08658 - Goat Anti-USP19 Antibody. It is currently untested in blocking **Product**

Comments applications.

Peptide

RSRNPGKDDCAKE Sequence

Ambient Instructions Shipping

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade