

EBP08705 - IMMUNIZING PEPTIDE

SKU: EBP08705



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_115972.2; NP_660187.1; NP_001357425.1
Immunogen	Peptide with sequence C-KVVEFDRPEVLRK, from the internal region (near C Terminus) of the protein sequence according to NP_115972.2; NP_660187.1; NP_001357425.1.
Product Comments	The peptide was used in the production of EB08705 - Goat Anti-MRP8 / ABCC11 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-KVVEFDRPEVLRK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

DOCUMENTS

- [Data Sheet](#)