

EBP08724 - IMMUNIZING PEPTIDE

SKU: EBP08724



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_004296.1; NP_647593.1; NP_647594.1; NP_647595.1; NP_647596.1; NP_647597.1; NP_647598.1; NP_647599.1; NP_647600.1; NP_647601.1
Immunogen	Peptide with sequence C-SDNAPAKGNKSPS, from the internal region of the protein sequence according to NP_004296.1; NP_647593.1; NP_647594.1; NP_647595.1; NP_647596.1; NP_647597.1; NP_647598.1; NP_647599.1; NP_647600.1; NP_647601.1.
Product Comments	The peptide was used in the production of EB08724 - Goat Anti-BIN1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-SDNAPAKGNKSPS
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)